

Dräger SSK 5000 Drug Testing

Ideal for surface drug testing applications, the Dräger SSK 5000 sampling system is quick and easy to use. When used with the Dräger DrugTest 5000, the sample can be transferred to a test cassette and screened for possible drug contamination. The sample can also be stored and transported for further testing.



D-62945-2009

Benefits

Surface Testing

The Dräger SSK 5000 consists of a vial containing a buffer solution and an integrated swab in the cap for sampling. The swab is used to wipe surfaces and test for drugs via the Dräger DrugTest® 5000 System.

Confirmation Testing

The substance of surfaces and unknown formulation which are transferred into the buffer solution using the SSK 5000 can be stored and transported safely in the sample container. For dosing of liquid, the sample container can be modified into a dropping bottle thanks to a special screw cap. The liquid to be tested can then be separated or reduced for subsequent analysis using droplet dosing.

Unknown Substances

Use the swab of the SSK 5000 to collect small amounts of, for example, an unknown granular formulation, transfer this into the buffer solution and test this for drugs via the Dräger DrugTest® 5000 System.

Related Products



D-54720-2012

Dräger DrugTest® 5000

No pipetting, no drips, no timing: Collecting an oral fluid sample is quick and easy with the Dräger DrugTest® 5000. The collected sample can then be analysed immediately for accurate results on the spot.

Related Products

ST-13226-2007



Dräger DCD 5000

Simple, fast and secure detection of intoxicants and drugs: The Dräger DCD 5000 does not only enable an easy collection of saliva, you can also store the samples and send them to the laboratory for analysis – without the sample being falsified along the way.

D-10639-2011



Dräger DrugCheck® 1200 STK

The reliable method to prevent drug abuse in a cost-effective way. Current drug influence can quickly and easily be detected with the saliva test.

Ordering Information

Dräger SSK 5000

Dräger SSK 5000, 20 Pieces	83 20 490
Pack sizes: 20 Pieces buffer-filled collection devices incl. swab-tap	

Dräger DrugTest 5000

Dräger DrugTest 5000 Analyzer, 1 Piece	83 19 900
Dräger DrugTest 5000 Test-Kits, 20 Pieces	83 19 970

CORPORATE HEADQUARTERS

Drägerwerk AG & Co. KGaA
 Moislinger Allee 53–55
 23558 Lübeck, Germany
www.draeger.com

Manufacturer:

Dräger Safety AG & Co. KGaA
 Revalstraße 1
 23560 Lübeck, Germany

REGION EUROPE CENTRAL AND EUROPE NORTH

Dräger Safety AG & Co. KGaA
 Revalstraße 1
 23560 Lübeck, Germany
 Tel +49 451 882 0
 Fax +49 451 882 2080
info@draeger.com

REGION EUROPE SOUTH

Dräger Safety France SAS
 25, rue Georges Besse
 92182 Antony Cedex
 Tel +33 1 46 11 56 00
 Fax +33 1 40 96 97 20
europe.south@draeger.com

REGION MIDDLE EAST, AFRICA

Dräger Safety AG & Co. KGaA
 Branch Office
 P.O. Box 505108
 Dubai, United Arab Emirates
 Tel +971 4 4294 600
 Fax +971 4 4294 699
contactuae@draeger.com

REGION ASIA PACIFIC

Draeger Safety Asia Pte Ltd
 25 International Business Park
 #04-20/21 German Centre
 Singapore 609916
 Tel +65 6308 9400
 Fax +65 6308 9401
asia.pacific@draeger.com

REGION CENTRAL AND SOUTH AMERICA

Dräger Panama S. de R.L.
 Complejo Business Park,
 V tower, 10th floor
 Panama City
 Tel +507 377 9100
 Fax +507 377 9130
contactcsa@draeger.com

Locate your Regional
 Sales Representative at:
www.draeger.com/contact

